

CLARIFLOW®-B

BIOLOGICAL RETENTIVE FILTER

Clariflow-B biological retentive filter cartridges are optimized to assure end-product purity and stability in the filtration of water, buffers, culture media, and biologicals used in the pharmaceutical and biopharmaceutical industry. Clariflow-B single layer of a mirrored-anisotropic polyethersulfone (PES) membrane produces a cartridge with high flow rates and consistent removal of microorganisms and inorganic particulates.

Clariflow-B is manufactured with polypropylene components and an inherently hydrophilic PES membrane that contains no added surfactants or wetting agents. The PES membrane also exhibits low protein-binding characteristics, is biologically inert, and is full bacterial retentive (LRV >7).

Each cartridge is fabricated in a cleanroom environment, pre-flushed with deionized water, and 100% integrity tested prior to shipment.



FEATURES & BENEFITS

- Bacterial retention validated in accordance with current version of ASTM F838 liquid bacterial challenge using *Brevundimonas diminuta*:
 - Economical sterilizing grade filtration membrane
- Single layer of Sterilizing 0.2um PES membrane design:
 - Sterilizing grade filtration validated with insert endcaps
- Integrity testing parameters correlated to bacterial retention:
 - For non-destructive integrity test to confirm sterilizing filtration
- Biologically inert and a low protein binding nature due to its construction materials:
 - Maximizes chemical compatibility and product yields
- PFAS Free construction:
 - Polypropylene components and PES filtration media

APPLICATIONS

Biological retention to sterilizing grade filtration of:

- | | |
|---------------------------|---------------------------------------|
| • Pharmaceutical Products | • Vaccines |
| • Parenteral Solutions | • Downstream Purification Bioproducts |
| • Ophthalmic Solutions | • Chromatography Column Buffers |
| • Diagnostic Reagents | • Blood and Plasma Fractionation |
| • Biological Growth Media | • Water for Injection |
| • Concentrated Proteins | |



TECHNICAL SPECIFICATIONS

Cartridge and Mini Cartridge:	
Filtration Membrane	Polyethersulphone
Upstream Support	Polypropylene
Downstream Support	Polypropylene
Inner Core	Polypropylene
Outer Cage	Polypropylene
End Caps	Polypropylene
Endcap Inserts	Polypropylene encapsulated 316L Stainless Steel, grafted Polysulphone
O-rings	Silicone, EPDM, Buna-N, Viton
Mini Capsules	
Body	Polypropylene
Connections	Polypropylene
Vent O-rings	Silicone, EPDM, Buna-N, Viton

EFFECTIVE FILTRATION AREA (EFA)

Size	Area (ft ²)	Area (m ²)
10 inch Cartridge	6.8	0.63
5 inch Cartridge	3.4	0.32
Mini Cartridge	2.5	0.23
Double Size Mini Capsule	2.5	0.23
Standard Size Mini Capsule	1.8	0.17
Half Size Mini Capsule	0.9	0.08


INTEGRITY TEST DATA

Rating	Bubble Point	Diffusional Flow	
μm	Min psi (bar)	Max value (ml/min)	Test Pressure psi (bar)
10" Cartridge			
0.1*	28* (1.9*)	40	55 (3.8)
0.2	45 (3.1)	20	32 (2.2)
0.45	30 (2.1)	20	24 (1.7)
0.6	17 (1.2)	20	13 (0.9)
0.8	15 (1.0)	20	11 (0.8)
5" Cartridge			
0.1*	28* (1.9*)	24	55 (3.8)
0.2	45 (3.1)	12	32 (2.2)
0.45	30 (2.1)	12	24 (1.7)
0.6	17 (1.2)	12	13 (0.9)
0.8	15 (1.0)	12	11 (0.8)
Mini Cartridge & Double Size Mini Capsule			
0.1*	28* (1.9*)	17	55 (3.8)
0.2	45 (3.1)	10	32 (2.2)
0.45	30 (2.1)	8	24 (1.6)
0.6	17 (1.2)	8	13 (0.9)
0.8	15 (1.0)	8	11 (0.8)
Standard Size Mini Capsule			
0.1*	28* (1.9*)	13	55 (3.8)
0.2	45 (3.1)	7	32 (2.2)
0.45	30 (2.1)	8	24 (1.6)
0.6	17 (1.2)	8	13 (0.9)
0.8	15 (1.0)	8	11 (0.8)
Half Size Mini Capsule			
0.1*	28* (1.9*)	6	55 (3.8)
0.2	45 (3.1)	4	32 (2.2)
0.45	30 (2.1)	8	24 (1.6)
0.6	17 (1.2)	8	13 (0.9)
0.8	15 (1.0)	8	11 (0.8)


*Values based on wetted membrane with IPA: H₂O (60%/40%) solution.

WATER FLOW RATES

10 inch Cartridge



Mini Cartridge and Mini Capsules



Note: Water flow rate for a Mini Cartridge is the same as for a Double Size Mini Capsule, the above graph applies for both presentations.

MAXIMUM OPERATING CONDITIONS

Differential Pressure

Cartridge: 80 psid (5.5 bar) @ 75°F (24°C)

Mini Capsule: 80 psid (5.5 bar) @ 75° F (24°C)

Steam Sterilization

Cartridge: Maximum 30 cycles of 30 minutes @ 257°F (125°C)

Mini Capsule: Max 10 Autoclave cycles of 30 mins @ 268°F (131° C)

Bacterial Retention

Clariflow-B filters have been validated as bacterial retentive filters as a standard option; and fully sterilizing grade filters when used with encapsulated insert end caps. They provide sterile filtrate when challenged with a liquid bacterial culture in accordance with current revision ASTM F838.

Biological Safety

All construction material in Clariflow-B conform with current biocompatibility testing per UPS <87>, <88> and <665> for moderate risk.

Product Assurance

Product is packaged and sealed in a protective polyethylene bag within a controlled manufacturing environment. Every filter device is integrity tested and flushed with DI water before shipment.

Packaging Options

- Non-sterile as a standard option
- Pre-sterilized mini capsule by autoclave

ORDERING INFORMATION

Cartridges: 25 — 0 — — — B

Insert Style		End Fitting		Nominal Length		Filter Rating		O-Ring	
Code	Description	Code	Description	Code	Length	Code	Micron	Code	Material
1	None (STD)	2	226/Flat	05	5" (130 mm)	001	0.1	0	Buna N
5	Encapsulated 316L SS	3	222/Flat	10	10" (250mm)	002	0.2	1	EPDM
6	Encapsulated Polysulfone	7	226/Fin	20	20" (500mm)	004	0.45	2	Silicone
		8	222/Fin	30	30" (750mm)	006	0.65	4	Viton®
				40	40" (1000mm)	008	0.8	N	None

Mini Cartridges: 25 — M D B M M — —

Filter Rating		O-Ring	
Code	Micron	Code	Material
001	0.1	0	Buna N
002	0.2	1	EPDM
004	0.45	2	Silicone
006	0.65	4	Viton®
008	0.8	N	None

Mini Capsules: 25 — C B — —

Size		Inlet/Outlet		Filter Rating		Vent O-Ring		Standard Package Options	
Code	Description	Code	Length	Code	Micron	Code	Material	Blank	Non-Sterile
H	Half	B	1/4" Hose Barb	001	0.1	0	Buna N	S	Pre-Sterilized
S	Standard	H	1/2" Hose Barb	002	0.2	1	EPDM		
D	Double	S	1 1/2" Sanitary Flange	004	0.45	2	Silicone		
		D	Walther 3/8 Male	006	0.65	4	Viton®		
		P	1/4" NPT (Male)	008	0.8	N	None		
		J	1/2" NPT (Male)						
		K	5/8" NPT (Female)						
		G	1/2" Swagelok®						

Parker Hannifin Corporation

Bioscience & Water Filtration Division

2340 Eastman Avenue

Oxnard, California 93030 USA

805 604 3400 | bwf.oxn.support@support.parker.com

www.parker.com/bioscience

© 2025 Parker Hannifin Corporation
Bioscience & Water Filtration Division

